# VanishPoint® II Syringe

Automated Retraction Technology

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Manufacturing Location  
Jumin Bio-Technologies Co., Ltd.  
1888 Huhang Road,  
Jianghai Economic Zone  
Fengxian, Shanghai  
P. R. China

PATENTS—WWW.RETRACTABLE.COM/PATENTS

© 2022. Retractable Technologies, Inc. All rights reserved.

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


Strzykawka VanishPoint® II

### Precautions:

- Single use only. Reuse of this device may result in exposure to bloodborne pathogens, including Hepatitis B virus (HBV), Hepatitis C virus (HCV), and human immunodeficiency virus (HIV).
- Contents are sterile, non-toxic, and non-pyrogenic. Do not use if product or package is damaged.
- Not made with natural rubber latex.
- Use only with attached needle. Needle cannot be changed.
- Automated needle retraction occurs only when barrel is emptied and plunger is fully depressed.
- For applications where full dose is not administered, expel remaining contents according to institutional policy and activate needle retraction.
- U.S. Federal Law restricts this device to sale by or on the order of a physician.

### Instructions for Use:

1. Prepare and give injection using aseptic technique according to institutional policy.
2. Continue depressing plunger to activate automatic needle retraction **while needle is still in patient. Full dose is administered when the plunger handle is fully depressed**, which activates the needle retraction mechanism.
3. Needle will automatically retract into syringe, preventing exposure to contaminated needle and rendering syringe non-reusable. In the event that needle retraction mechanism does not activate, discard syringe in an appropriate sharps container per protocol of institution. Do not recap contaminated needles.
4. Dispose of VanishPoint® II syringe in an appropriate sharps container per protocol of institution.

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EU Auth. Rep. (MDD93/42/EEC) MDSS GmbH Schiffgraben 41 D-30175 Hannover Germany