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







VanishPoint® Blood Collection System

VanishPoint® Blutentnahmesystem
Système de prélèvement sanguin VanishPoint®
Sistema per raccolta di sangue VanishPoint®
Sistema de extracción de sangre VanishPoint®
VanishPoint® Bloedverzamelingsysteem

Sistema de coleta de sangue VanishPoint®
VanishPoint® 採血システム
VanishPoint® Kan Alma Sistemi
System pobierania krwi VanishPoint®

VanishPoint® Blood Collection Tube Holder
VanishPoint® Blutentnahmeröhrchenhalter
Porte-tubes de recueil de sang VanishPoint®
Porta-fiale per raccolta di sangue VanishPoint®
Porta-tubo para muestras sanguíneas VanishPoint®
VanishPoint® Bloedverzamelingsbuis houder
Adaptador para tubos de coleta de sangue VanishPoint®
VanishPoint® 採血チューブホルダー
VanishPoint® Kan Alma Tüp Tutucusu
Uchwył na próbówki do pobierania krwi VanishPoint®

VanishPoint® Small Diameter Tube Adapter
VanishPoint® Adapter für Röhrchen mit kleinem Durchmesser
Adaptateur pour tubes de petit diamètre VanishPoint®
Adattatore per fiale di piccolo diametro VanishPoint®
Adaptador para tubos de diámetro pequeño VanishPoint®
VanishPoint® Adapter voor kleine buisdiameter
Adaptador para tubos de pequeno diámetro VanishPoint®
VanishPoint® 小径チューブアダプター
VanishPoint® Küçük Çaplı Tüp Adaptörü
Adapter próbówki o malej średnicy VanishPoint®

REAL SAFETY™

Corporate Headquarters
Retractable Technologies, Inc.
511 Lobo Lane, P.O. Box 9
Little Elm, Texas 75068-0009
www.retractable.com

Sales and Marketing
Retractable Technologies, Inc.
511 Lobo Lane, P.O. Box 9
Little Elm, Texas 75068-0009
Fax: (972) 294-4400
Toll Free 1-888-703-1010

Manufacturing Location
Retractable Technologies, Inc.
511 Lobo Lane, P.O. Box 9
Little Elm, Texas 75068-0009

EC REP EU Auth. Rep.
MDSS GmbH
Schiffgraben 41
D-30175 Hannover
Germany

PATENTS--WWW.RETRACTABLE.COM / PATENTS

© 1999, 2023. Retractable Technologies, Inc. All rights reserved.

Precautions:

- Single use only. Reuse of this device may result in exposure to bloodborne pathogens, including Hepatitis B virus (HBV), Hepatitis C virus (HCV), and human immunodeficiency virus (HIV).
- Do not use if product is damaged.
- This device is not made with natural rubber latex.
- For use with conventional multiple sample blood collection needles up to 1 1/2" (38 mm) in length.
- Not for use with winged needle sets or luer adapters.
- Do not attempt to remove blood collection needle from tube holder after retraction.
- U.S. Federal Law restricts this device to sale by or on the order of a physician.

Instructions for Use:

1. Securely thread multiple sample blood collection needle into tube holder. Remove needle cover when ready to use.
2. Perform venipuncture per protocol of institution.
3. Insert blood collection tube into the tube holder, pushing forward to start blood flow. Remove blood collection tube when filled and repeat with additional tubes as required. For small diameter tubes, use with VanishPoint® small diameter tube adapter.
4. Once last tube has been removed and **needle is still in patient's vein**, close end cap securely into tube holder.
5. Needle will automatically retract from patient's vein into tube holder. In the event that needle retraction mechanism does not activate, discard tube holder in an appropriate sharps container per protocol of institution. Do not recap contaminated needles.
6. Dispose of tube holder in an appropriate sharps container per protocol of institution.

Precautions:

- Single use only. Reuse of this device may result in exposure to bloodborne pathogens, including Hepatitis B virus (HBV), Hepatitis C virus (HCV), and human immunodeficiency virus (HIV).
- This device is not made with natural rubber latex.
- For use with small diameter blood collection tubes and VanishPoint® blood collection tube holder only. See tube holder instructions for use.
- U.S. Federal Law restricts this device to sale by or on the order of a physician.

Instructions for Use:

1. Securely thread multiple sample blood collection needle into tube holder.
2. Perform venipuncture per protocol of institution.
3. Insert small diameter blood collection tube into tube adapter.
4. Insert tube adapter into tube holder, pushing forward to start blood flow.
5. When blood collection tube is filled, remove tube adapter from tube holder and blood collection tube from tube adapter. Repeat with additional tubes as required.
6. Dispose of tube adapter per protocol of institution.

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