



U-100 INSULIN SAFETY SYRINGES

Description	Needle Size	# of Units per Case	# of Boxes per Case	Catalog Number
0.5mL VanishPoint U-100 Insulin Syringe (half-unit)	30G x 3/16" (5mm)	800	8	15391
0.5mL VanishPoint U-100 Insulin Syringe (half-unit)	30G x 5/16" (8mm)	800	8	15371
0.5mL VanishPoint U-100 Insulin Syringe	30G x 3/16" (5mm)	800	8	15291
0.5mL VanishPoint U-100 Insulin Syringe	30G x 5/16" (8mm)	800	8	15271
0.5mL VanishPoint U-100 Insulin Syringe	30G x 1/2" (12.7mm)	800	8	15221
1mL VanishPoint U-100 Insulin Syringe	30G x 3/16" (5mm)	800	8	10291
1mL VanishPoint U-100 Insulin Syringe	30G x 5/16" (8mm)	800	8	10271
1mL VanishPoint U-100 Insulin Syringe	29G x 5/16" (8mm)	800	8	10281
1mL VanishPoint U-100 Insulin Syringe	29G x 1/2" (12.7mm)	800	8	10211

## ADVANTAGES

- Accurate dose delivery
- Designed for patient comfort
- Easy to use
- Safe and effective
- Both hands remain behind the sharp
- Reduces risk of needlestick injury
- Non-reusable when activated
- Remains safe through disposal
- Minimizes disposal volume

### Other available safety products



EasyPoint®  
Needle

VanishPoint®  
Blood Collection  
Set

EasyPoint®  
Blood Collection  
Tube Holder  
with Needle

VanishPoint®  
IV Catheter



**RETRACTABLE**  
TECHNOLOGIES, INC.

[www.retractable.com](http://www.retractable.com)

511 Lobo Lane • Little Elm • Texas • 75068-5295 • USA  
Tel.: (972) 294-1770 • Fax: (972) 294-4400 • Toll Free: 1-888-703-1010

©2023 Copyright Retractable Technologies, Inc. All rights reserved. M161 REV. 0 EN, JAN 2023



# Safe and accurate



**RETRACTABLE**  
TECHNOLOGIES, INC.



**VanishPoint**<sup>®</sup>  
U-100 INSULIN SAFETY SYRINGES

**Real  
Safety**<sup>™</sup>

## NEEDLESTICK PREVENTION

- Safe and effective
- Reduces risk for needlestick injury
- Both hands remain behind the sharp



Automated pre-removal  
retraction effectively  
reduces the risk  
of needlestick  
injuries.

Integrated safety mechanism allows for optimum view of syringe scale, promoting dose preparation ease & accuracy

VanishPoint<sup>®</sup> Insulin syringes are available with ½ unit increments for precise small dose delivery

VanishPoint<sup>®</sup> Insulin syringes are available with 5mm, 8mm and 12.7mm needles to ensure subcutaneous insulin delivery and maximize patient comfort