


Before After

Available with
5mm needles!



0.5mL VanishPoint[®]
U-100 Insulin Syringe

SAFE AND ACCURATE



Specifications

- Automated pre-removal retraction effectively reduces the risk of needlestick injuries
- Integrated safety mechanism allows for optimum view of syringe scale, promoting dose preparation ease & accuracy
- VanishPoint[®] Insulin syringes are available with 5mm, 8mm and 12.7mm needles to ensure subcutaneous insulin delivery and maximize patient comfort
- VanishPoint[®] Insulin syringes are available with 1/2 unit increments for precise small-dose delivery

Advantages

- Accurate dose delivery
- Designed for patient comfort
- Easy to use
- Safe and effective
- Both hands remain behind the sharp
- Reduces risk of needlestick injury
- Non-reusable when activated
- Remains safe through disposal
- Minimizes disposal volume



RETRACTABLE[®]
TECHNOLOGIES, INC.

511 Lobo Lane • Little Elm • Texas • 75068-5295 • USA
Tel: (972) 294-1770 • Fax: (972) 294-4400
Toll Free: 1-888-703-1010



www.retractable.com